



Manufacturer & Importer of Dental Materials & Equipment

International Catalog 2026



MFG/MD/2023/000109
MFG/MD/2023/000100



Prime-Dental.in



+91 9619635313



info@prime-dental.in

About us



Company Overview:

Founded in 1987 by Mr. Bharat Thakkar and Dr. Paresh Thakkar, Prime Dental Products Pvt. Ltd. has grown from a dental trading enterprise into one of India's trusted manufacturers and importers of premium dental medical devices.

Prime Dental operates from a world-class, ISO 13485:2016 certified manufacturing facility spread across 13,000 sq. ft. in MIDC Wagle Estate, one of Western India's leading industrial hubs.

With a strong reputation for quality, innovation, and reliability, we serve the dental fraternity across India through an extensive dealer network. Our range of dental products is also exported internationally through valued global trade partners.

Committed to advancing dentistry through superior products and long-term partnerships, Prime Dental continues to strengthen its legacy under the leadership of the next generation, Mr. Tanay Thakkar and Mr. Jay Thakkar, who have been actively driving the company's growth and vision for over a decade.

Dental Materials - 'Make in India' initiative:

We are committed to serve the dental community offering globally acceptable quality products at competitive prices and timely delivery.

Our product range include;

- **Endodontic materials:** MTA Bioceramic Sealer, RC Help (EDTA Gel), RC Cal (Calcium hydroxide in ready to use paste form for temporary RC filling), RC Seal (Epoxy Bonded RC Sealant).
- **Prosthodontic materials:** Prime Chrome (Chromatic alginate impression material), Retra Gel (Retraction Gel) & the foremost manufacturer of Addition silicone impression material (Ad Sil Acura), Bite Registration Putty.
- **Conservative materials:** Restorite Etching Gel (37% orthophosphoric acid gel), TMP-RS (Ready to use temporary filling putty material), Dent Wash (Liquid EDTA).
- **Finishing materials:** Aster prophyl (pumice based prophylaxis paste) and Aster Compo (Alumina based composite polishing paste).
- **Composites:** Restorite range of Composites offers first in class physical properties viz. Bulk Fill, Flow Viscous, Pit & Fissure Sealant, Micro Hybrid Composite, Nano Hybrid, Core Build-up, Restorite Bond 5G (Total etch), Restorite Bond 7G (Self etch).

Import & Distribution:

Since 1991, Prime Dental is representing many renowned international brands in India.

Dental Equipment:

Since 2007, We are dealing in the finest quality of fully programmable electrical Dental chairs & heavy/light Suctions. We have various models of diverse range to cater all your requirements. We specialize in new clinic setup & up-gradation with robust after-sale service.

OBL Products:

Prime Products: Prime X-Ray Film, Prime Cord (Knitted Retraction cord), Prime Disposables (Gloves, Face Mask, Saliva Ejectors).

Allure: Airtors.

Achievements:

Best Dental Manufacturer - Make in India
at 1st National Oral Healthcare Sushruta Awards
hosted by Indian Dental Association (IDA) Foundation
at the hands of Honourable Former President of India,
Dr. A. P. J. Abdul Kalam



Manufactured Products

Page

For Irrigation 01

Prime Chlorhexidine 2% 500ml **2026 Launch**

RC Twents 10Pcs (30 Gauge)

EDTA 01

RC Help 3g / 2 x 3g / 20g

Dent Wash 30ml

Related Product

Citra Wash 30ml

Calcium Hydroxide 02

BioTemp 3g

RC Cal 2g / 2 x 2g

RC Cal-I 2g / 2 x 3g (with Iodoform)

RC Pex 3g (with Iodoform)

Restorite Cal 2g (Light Cured & Rigid Setting)

Temporary Filling Calcium Hydroxide Powder

Temporary Filling Material 02

Temporary Filling Eugenol Liquid 15ml

Temporary Filling Eugenol Liquid 110ml

Temporary Filling Zinc Oxide Powder 100g

TMP-RS 40g Eugenol Free

Zinc Oxide Eugenol Cement Quick Set

Polymer Reinforce (IRM) 15g + 10ml

Root Canal Sealers 03

Prime MTA Bioceramic Root Canal Sealer 3g

RC Fill 15g Powd + 10ml Liq + 5g Acc

RC Fill-I 15g Powd + 10ml Liq + 5g Acc (with Iodoform)

RC Seal 10g Powd & 10ml Liq

RC Seal Plus (12g Liq + 3.5 Liq)

Root Repair Material 03

Prime MTA BioPutty 2g

Prime MTA Packable 1g Powd + 5ml Liq

Page

Solvents of Obturation Material 04

RC Soft 15ml

RC Solve 15ml

Root Canal Disinfectant 04

Primecresol 15ml

Prime CMCP 15ml **2026 Launch**

Prime DeVitex (3g) - Devitalizing Paste **2026 Launch**

Root Surface Planning 04

Root Help Gel (3g) - 24% EDTA **2026 Launch**

Root Help Liquid (15ml) - 15% Citric Acid **2026 Launch**

Hemostat & Tissue Management 05

Dent Retract 30ml

Retra Gel 3ml (Hemostat Gel)

Retra Liquid 30ml (Hemostat Liquid)

A-Silicone (Ad Sil Acura) 06

Bite Registration Putty 2 x 100ml

Putty - Hard Set 800g (Base + Catalyst)

Putty - Soft Set 800g (Base + Catalyst)

Light Cartridge 2 x 50ml (with Tips)

Light Tubes 2 x 25ml

Alginate 05

Algingum 500g

Prime Chrome 450g

Related Product

Clean Tray 250g

Impression Material 06-07

Green Stick 10 Sticks (Approx 90g)

Image 125g Base & 85g Catalyst

Impression Composition 5 Cakes (Approx 200g)

Temporary Luting 07

Templute 50g Base & 15g Catalyst

Manufactured Products

Page

Page

Glass Ionomer Cement (GIC) 07

Prime Cem Ceramic Reinforced

15g Powd+ 10ml Liq GIC

Prime Cem Luting Type 1

15g Powd+ 10ml Liq GIC

Prime Cem Restorative Type 2

15g Powd+ 10ml Liq GIC

Restorite Ionomer Flow (2g) - LC GIC **2026 Launch**

Etchant 08

Etching Gel 3ml / 2 x 3ml / 2 x 25ml

(Blue colour - Flowable)

Restorite Etching Gel 3ml / 2 x 3ml / 2 x 25ml

(Purple colour - Thixotropic)

Etching Liquid 15ml

Bonding Agent 08

Restorite Bond 5G 5ml

Restorite Bond 7G 5ml

Flowable Composite 09

PurpleSeal - Purple Colour Composite 2g

Restorite Bulk Fill (Flowable) 2g

Restorite Bulk Fill Kit (Flowable) 2g x 4

Restorite Flow Viscous 2g

Restorite Flow Viscous Kit 2g x 4

Restorite Ionomer Flow (2g) - LC GIC **2026 Launch**

Restorite Pit & Fissure Sealant 2g

Packable Composite 10

BiteRise (4g) - Blue Colour Composite **2026 Launch**

Restorite Bulk Fill (Packable) 4g

Restorite Bulk Fill Kit (Packable) 4g x 7

Restorite Core Build-up Translucent 4g (6mm depth of cure)

Restorite Micro Hybrid 4g

Restorite Micro Hybrid Kit 4g x 7

Restorite Nano Hybrid 4g

Restorite Nano Hybrid Kit 4g x 7

Liquid Resin 11

Restorite Cavo Surface Sealer 5ml

Restorite Modelling Resin 5ml

Orthodontic Material 11

BiteRise (4g) - Blue Colour Composite **2026 Launch**

Restorite OrthoGrip Composite 5g

Restorite OrthoGrip Primer 5ml

Restorite OrthoGrip Composite Kit

Gingival Barrier 12

GingiSafe 1.2ml

GingiSafe Kit 1.2ml x 4

Bleaching System 12

LilyBright **2026 Launch**

Desensitisers 12

Desensitiser 3g / 2 x 3g

Restorite Desensitiser 5ml

Prophylaxis & Polishing System 12-13

Aster Compo 4g

Aster Prophy 45g

(Lemon/Orange/Strawberry/Mint)

Aster Prophy 4 x 45g (Assorted Flavours)

Miscellaneous 13

BioFilmX 500ml

Caries Indicator 15ml

Debubblizer 100ml Spray / 500ml Refill Jar

Dent Dry 30ml

NT Fog 100ml

ProstoClean 5ml **2026 Launch**

Twin Packs available in



For Irrigation

Prime Chlorhexidine 2026 Launch



2%


2% Chlorhexidine Digluconate Solution for Effective Irrigation during Endodontic Procedures



I Open cap Remove aluminium seal



II Puncture syringe adaptor with needle on three sides (Do not puncture centre part first)



III Puncture syringe adaptor from centre hole with needle (Straight puncture is recommended)



IV Engage the syringe firmly on syringe adaptor and start pulling syringe piston

Available Packing : 500ml Bottle in 2%

RC Twents



30 Now in gauge

Irrigation Needle

- Irrigation needle with two side vents.
- Side vent avoids irrigation fluid to be pushed apically.

Available Packing : 10pcs in a pack

EDTA

RC Help



Lubricating Paste for RC Biomech.

- EDTA paste.
- Contains carbamide peroxide.
- Increases RCT instrument's life.
- Gentle O₂ release while irrigating with NaOCl.
- Cleanses canal

Available Packing • **Regular Pack** : 3g syringe & 5 application needles
 • **Twin Pack** : 2 x 3g syringe & 10 application needles
 • **Big Pack** : 20g big syringe, 25 application needles, 1 contactor & 3g empty syringe

Dent Wash



Aqueous EDTA

- 17% EDTA
- For smear layer removal.
- Use before bonding of dentin.
- Use before root canal sealant.

Available Packing : 30ml Bottle

Related Product

Citra Wash



20% Citric Acid Solution

- Cleansing/conditioning of prepared canals, and
- Dissolving/cleaning residues of calcium hydroxide
- Removing smear layer

Available Packing : 30ml Bottle

Calcium Hydroxide

BioTemp



Easy & Ready-to-Use Bioceramic Intra Canal Temporary Filling Material

- Radiopaque
- Intracanal dressing material
- Non-setting, ready-to-use paste
- Contains silicates for gradual release of calcium ions
- Can be left in the root canal for up to 3 weeks
- Easy to remove
- Recommended for use in cases of perforation & resorption prior to Prime's MTA Putty, TA Packable & MTA Bioceramic Sealer

Available Packing : 3g Syringe & 5 Capillary Tips

RC Cal-I



With Iodoform

Ready to Use Paste. Intracanal Temporary Filling

- Contains barium sulphate. Radio opaque.
- Highly alkaline Ca(OH)₂ paste.
- Can be easily removed by reamers.
- Non drying paste.

Available Packing • **Regular Pack** : 2g syringe & 5 application needles
 • **Twin Pack** : 2 x 2g syringe & 10 application needles

Restorite Cal



Light cured Calcium Hydroxide Cavity Base/Liner

- Easy to use, light curable, one component system
- Used for direct and indirect pulp capping
- Used as cavity liner under filling material
- Low solubility

Available Packing : 2g syringe

RC Cal



Without Iodoform

Ready to Use Paste. Intracanal Temporary Filling

- Contains barium sulphate. Radio opaque.
- Highly alkaline Ca(OH)₂ paste.
- Can be easily removed by reamers.
- Non drying paste.

Available Packing • **Regular Pack** : 2g syringe & 5 application needles
 • **Twin Pack** : 2 x 2g syringe & 10 application needles

RC Pex




Calcium Hydroxide Iodoform Paste Temporary Root Canal Filling Material

- Highly radiopaque calcium hydroxide paste with iodoform.
- Capillary tips for easy dispensing.

Available Packing : 3g syringe & 5 Capillary Tips

Temporary Filling Calcium Hydroxide Powder



Powder For Dental Cement For Cavity Base

- It can be used for pulp capping. Need to mix with water for injection to make paste of desired consistency while using

Available Packing : 10g Bottle

Temporary Filling Material

Temporary Filling Eugenol Liquid



- Used with zinc oxide powder for temporary filling.

Available Packing : 15ml Bottle / 110ml Bottle

TMP-RS



Temporary Filling Material

- Eugenol Free
- Temporary filling for cavities.
- Sets with moisture.
- Easy to place and remove.
- Two wads to offer maximum protection from moisture.

Available Packing : 40g Bottle

Temporary Filling Zinc Oxide Powder



- Used with eugenol liquid for temporary filling material.

Available Packing : 100g Bottle


Zinc Oxide Eugenol Cement Quick Set Polymer Reinforced (Intermediate Restorative Material)



Packing : 15g Powder / 10ml Liquid

Root Canal Sealers

Prime MTA Bioceramic Root Canal Sealer



- Ready-to-Use Paste.
- Premixed, Easy-to-Use.
- Radiopaque
- Bioactive & Biocompatible
- Hydrophilic

Available Packing : 3g Paste & 15 Capillary Tips



RC Fill
Without Iodoform

Root Canal Sealant & Filling
 Contains Calcium Hydroxide, Zinc oxide & Eugenol.
 Radio opaque.
 Bioresorbable set mass.

Available Packing : 15g Powder, 10ml Liquid & 5g accelerator



Available Packing : 15g Powder, 10ml Liquid & 5g accelerator

RC Seal

Epoxy Bonded Root Canal Sealant



- Root canal sealant.
- Epoxy bonded.
- Virtually formaldehyde free.
- Shrink free set.
- Eugenol free.
- Shall not retard composite setting.
- Non Bio Resorbable.

Available Packing : 10g Powder & 10ml Liquid
 ● 10ml refill liquid also available separately

RC Seal Plus

Epoxy Bonded Root Canal Sealant




- Hexamine free Epoxy bonded root canal sealant.
- Easy to mix paste-paste system in a two-barrel syringe.
- High radio opacity.
- Long term sealing properties.
- Good dimensional stability.

Available Packing : 1 x 15.5g Dual syringe

Root Repair Material

Prime MTA BioPutty

Ready-to-Use Bioceramic Putty for Root Repair



- Putty Consistency
- Non-Staining
- Highly Radiopaque

Available Packing : 2g Syringe

Prime MTA Packable

MTA Packable can be used for Apexification, Repairing root perforations, Internal root resorption and Pulp capping.



Available Packing : 1g Powder & 5ml Liquid

Solvents of Obturation Material

RC Soft

Softens Resin Based Cements



- Softening agent for epoxy resin based root canal sealants.

Available Packing : 15ml bottle

RC Solve

GP & ZnOE Solvent



- Softens GP points.
- Softens zinc oxide cement.
- Controlled penetration to avoid disturbing apical seal.
- Shall not retard composite setting.

Available Packing : 15ml bottle

Root Canal Disinfectant

Primecresol



Formocresol

Available Packing : 15ml Bottle

Prime CMCP

2026 Launch



Camphorated Parachlorophenol

- Disinfectant Antimicrobial Agent for Root Canal
- Corticosteroid-Free
- Antiseptic Action by Camphor
- Bactericidal Action by Para-Chlorophenol
- For Professional Use Only

Available Packing : 15ml Bottle

Prime DeVitex

2026 Launch



Pulp Devitalising Paste. Arsenic Free.

- Bright Red to Easily Identify
- For Pulp Necrosis & Devitalization
- Fast-Acting

Available Packing • 3g syringe

Root Surface Planning

Root Help Gel

2026 Launch

24% EDTA
For Bio Root Planning

- Chelating Agent during Root Planning
- Removes Smear Layer, Exposes Collagen for Improved Periodontal Ligament Cell Adhesion



Available Packing • 3g syringe

Root Help Liquid

2026 Launch

15% Citric Acid
Root Planning Liquid

- Disinfectant Antimicrobial Agent for Root Canal
- Corticosteroid-Free
- Antiseptic Action by Camphor
- Bactericidal Action by Para-Chlorophenol
- For Professional Use Only



Available Packing : 15ml Bottle

Hemostat & Tissue Management

Dent Retract



Aqueous Ferric Chloride

- Helps in retraction of gingival.
- Does not retard the rubber base.

Available Packing : 30ml Bottle

Retra Gel

25% Aluminum Chloride Gel



- Retra Gel is a 25% Aluminum Chloride aqueous viscous gel.
- No residue or stain making it ideal for use in the esthetic zone.
- The solutions facilitate the insertion of the cord into the sulcus



Available Packing : 3ml Syringe & 5 application needles

Retra Liquid


25% Aluminium Chloride Aqueous Astringent Solution



- For chemical retraction of gingival margin
- Clear transparent solution
- Non-staining
- Use along with Prime Retraction Cord
- Wash after use before final impression

Available Packing : 30ml Bottle

Alginate



Algingum CE
 Non Chromatic



Available Packing : 500g pouch

REF : ALG001

Prime Chrome CE
 Chromatic

Thixotropic
 Rubber effect
 Dust free
 Easy mix formula
 Stable for 50 hours under humidity



Available Packing : 450g Pouch

REF : PCR001



Related Product

Clean Tray CE



Alginate and Plaster Residue Cleaner

- 250g powder for 2.5 litre solution.
- Biodegradable powder.
- Can clean up to 125 trays.

Available Packing : 250g Powder

REF : CTA001

A-Silicone (Ad Sil Acura)

Ad-Sil Acura Bite Registration Putty **CE**



- No Heating required
- Addition Silicone based Bite Registration Putty
- Shore A - 87
- Easy-to-Use
- Can be used for indexing
- No wastage compared to Bite Registration in Cartridge form

Available Packing : 100ml Base & 100ml Catalyst

REF : ABR001

Hard Set **CE**



Shore A - 80



- Addition silicone impression material.
- Accurate detailed reproduction.
- Highly stable impression.
- Compatible to all dye materials.
- Suitable for all impression technique.
- Ideal choice of impression for multiple pouring.

Available Packing : 400g Base & 400g Catalyst

REF : ASH001

Ad-Sil Acura Light Cartridge **CE**



Polyvinyl Siloxane Impression Material

- Shore A-65.
- Light body available in dual chamber cartridge.
- Hydrophilic nature.
- Can be used with any other addition silicone putty.

Available Packing : Cartridge 2 x 50ml (Base & Catalyst) with Tips

REF : ALC001

Ad-Sil Acura Putty **CE**

Soft Set **CE**



Shore A - 60

Available Packing : 400g Base & 400g Catalyst

REF : ASP001

Ad-Sil Acura Light Tubes **CE**



Polyvinyl siloxane impression material

- Shore A-65.
- Addition silicone impression material.
- Accurate detailed reproduction.
- Highly stable impression.
- Compatible to all dye materials.
- Suitable for all impression technique.
- Ideal choice of impression for multiple pouring.

Available Packing : Tubes 25ml Base & 25ml Catalyst

REF : ALT001

Impression Material

Image **CE**



Impression Paste for Edentulous Cases

- Zinc oxide base.
- Eugenol free.
- No irritation to lip & mucosa.

Available Packing : 125g Base & 85g Catalyst

REF : IMG001

Green Stick CE

Tracing Sticks



- Low fusing compound.
- Soft green tracing sticks.
- For border moulding work.

Available Packing : 10 Sticks in a pack (90g)

REF : GRS001

Impression Composition CE

For Accurate Edentulous Impressions



- High fusing compound.
- Functional impression compound.
- Accurate impression.
- Suitable for primary impression of complete denture prosthesis.

Available Packing : 5 cakes in a pack (200g)

REF : IPC001

Temporary Luting

Templute

Temporary Luting for C & B



- Two paste system.
- Eugenol free.
- Shall not retard any bonding agent.
- C & B can be removed easily.

Available Packing : 50g Base & 15g Catalyst

Glass Ionomer Cement (GIC)

Prime Cem Luting

Glass Ionomer Cement



Available Packing : 15g Powder, 10ml Liquid, Scoop/Spatula & Mixing Pad

Prime Cem Restorative

Glass Ionomer Material



Available Packing : 15g Powder, 10ml Liquid, Scoop/Spatula & Mixing Pad

Prime Cem Ceramic Reinforced

Glass Ionomer Material
 Contains Zirconium Oxide



Available Packing : 15g Powder, 10ml Liquid, Scoop/Spatula & Mixing Pad

Restorite Ionomer Flow

2026 Launch

Light Cure Glass Ionomer
 Composite Cement



- High Compressive Strength
- Radiopaque
- Composition: Acrylate Monomer Glass Ionomer-Based Filler



Available Packing : 2g syringe

Etchant

Etching Gel

Blue Color

Flowable



37% Phosphoric Acid
 Enamel Etchant.
 Silica Free

Available Packing • **Regular Pack:** 3ml syringe & 5 application needles • **Twin Pack:** 2 x 3ml syringe & 10 application needles
 • **Big Pack:** 2 x 25ml syringe, 50 application needles, 1 contactor & 2 x 3ml empty syringe.

Restorite Etching Gel

Purple Color

Thixotropic



37% Phosphoric Acid
 Enamel Etchant.
 Contains Silica

Available Packing • **Regular Pack:** 3ml syringe & 5 application needles • **Twin Pack:** 2 x 3ml syringe & 10 application needles
 • **Big Pack:** 2 x 25ml syringe, 50 application needles, 1 contactor & 2 x 3ml empty syringe.

Etching Liquid



37% Orthophosphoric acid

- Enamel Etchant.
- Silica free.
- Thick liquid consistency.

Available Packing : 15ml Bottle

Bonding Agent

Restorite Bond 5G

5th Generation Bonding Agent for
 Dentine /Enamel



- Light-cured one component Total Etch System.
- Provides very good adhesion to dentine and enamel.

Available Packing : 5ml bottle

Restorite Bond 7G

7th Generation Self Etch
 Bonding Agent For Dentine & Enamel



- Light cured one component Self Etch Bonding System
- Provides very good adhesion to dentine and enamel

Available Packing : 5ml bottle

Flowable Composite

PurpleSeal



Flowable Composite for Sealing Orifice and Pulpal Floor

- Light-Cured Flowable Purple Composite
- High Filler Load: 67%
- Easy Identification of Orifice During Retreatment
- Can Also Be Used for Pulpal Floor Sealing
- Radiopaque

Available Packing : 2g Syringe & 5 Applicator Tips

Restorite Bulk Fill (Flowable)

Light Cure Flowable Bulk Fill Composite Filling Material



- Depth of cure 4mm
- Available Shade
Universal



Available Packing : 2g syringe

Restorite Flow Viscous

Light Cure Flowable Viscous Composite Filling Material



Available Shade
A1
A2
A3
B2



Available Packing : 2g syringe

Restorite Ionomer Flow

2026 Launch

Light Cure Glass Ionomer Composite Cement

- High Compressive Strength
- Radiopaque
- Composition: Acrylate Monomer Glass Ionomer-Based Filler



Available Packing : 2g syringe

Restorite Bulk Fill Kit (Flowable)



Contents:

- Restorite Bulk Fill Flowable (2gm) - 4
- Application Needles for Restorite Bulk Fill Flowable - 20



Shades included : Universal - 4

Restorite Flow Viscous Kit



Contents:

- Restorite Flow Viscous (2gm) - 4
- Application Needles for Restorite Flow Viscous - 20



Shades included: A1, A2, A3 & B2

Restorite Pit & Fissure Sealant

Light Cured Universal Shade Fissure Sealant

Available Shade
Universal Shade



Available Packing : 2g syringe

Packable Composite

BiteRise

2026 Launch



Bite Raising Composite

- for Easy Identification
- For Creating Temporarily Bonded, Raised Occlusal Surfaces of Disoccluded Teeth
- Can Also Be Used for Deep bite and Cross bite
- Easy to Apply and Remove



Available Packing : 4g syringe

Restorite Bulk Fill (Packable)

Light Cure Packable Bulk Fill Composite Filling Material

• Depth of cure 4mm

Available Shade
Universal



Available Packing : 4g syringe

Restorite Micro Hybrid

Light Cure Resin-Based Dental Restorative Material

Available Shade
A1, A2, A3, A3.5, B1 & B2



Available Packing : 4g syringe

Restorite Nano Hybrid

Light Cure Resin-Based Dental Restorative Material

Available Shade
A1, A2, A3, B2



Shades included : 4g Syringe

Restorite Core Build-Up

• Depth of Cure 6mm



Available Shade
Translucent Shade



Available Packing : 4g syringe

Restorite Bulk Fill Kit (Packable)



Contents:

- Restorite Bulk Fill Packable (4gm) - 7
- Restorite Bond 7G (5ml) - 1
- Restorite Etching Gel (3ml) - 1
- Aster Compo (4g) - 1
- Application Needles for Restorite Etching Gel - 5
- Microbrush - 10

Shades included: Universal - 7

Restorite Micro Hybrid Kit



Contents:

- Restorite Micro Hybrid (4gm) - 7
- Restorite Bond 7G (5ml) - 1
- Restorite Etching Gel (3ml) - 1
- Aster Compo (4g) - 1
- Application Needles for Restorite Etching Gel - 5
- Microbrush - 10

Shades included: (A1 - 2) (A2 - 2) (A3 - 2) & (B2 - 1)

Restorite Nano Hybrid Kit



Contents:

- Restorite Nano Hybrid (4gm) - 7
- Restorite Bond 7G (5ml) - 1
- Restorite Etching Gel (3ml) - 1
- Aster Compo (4g) - 1
- Application Needles for Restorite Etching Gel - 5
- Microbrush - 10

Shades included: (A1 - 2) (A2 - 2) (A3 - 2) & (B2-1)

Liquid Resin

Restorite Cavo Surface Sealer



Light Cured; Unfilled Methacrylate Resin-Based Cavo Surface Sealer

- For better margin integrity.
- For filling micro cracks.
- Can be used for GI and Composites.

Available Packing : 5ml bottle

Restorite Modelling Resin



Composite Sculpting Agent

- Low viscosity light curable resin for sculpting of composite.

Available Packing : 5ml bottle

Orthodontic Material

BiteRise
 2026 Launch



Bite Raising Composite

- **Blue Color** for Easy Identification
- For Creating Temporarily Bonded, Raised Occlusal Surfaces of Disoccluded Teeth
- Can Also Be Used for Deep bite and Cross bite
- Easy to Apply and Remove

Available Packing : 4g syringe

Restorite OrthoGrip Composite

2026 Launch

Light Cure Orthodontic Bracket Adhesive



Available Packing : 5g syringe

Restorite OrthoGrip Primer

2026 Launch



Available Packing : 5ml bottle

Restorite OrthoGrip Composite Kit

Light Cure Orthodontic Adhesive System



Available Packing • Kit

- Restorite OrthoGrip Composite (5g) - 1
- Restorite OrthoGrip Primer (5ml) - 1
- Restorite Etching Gel (3ml) - 2
- Application Needles for Restorite Etching Gel - 10
- Microbrush - 10



Gingival Barrier

GingiSafe

Light Cure Gingival Barrier / Liquid Dam



Available Packing : 1.2ml syringe

GingiSafe Kit

Light Cure Gingival Barrier / Liquid Dam



Available Packing : Kit : GingiSafe (1.2ml) - 4, Application Needles - 8

Bleaching System

LilyBright

2026 Launch



Fast Acting In-Office Whitening System 35% Hydrogen Peroxide with pH-Activation and Desensitiser

Available Packing • Kit

- LilyBright Powder (0.3g) - 2
- LilyBright Liquid (3ml) - 2
- GingiSafe (1.2ml) - 1
- Microbrush - 2
- Applicator Tips - 2



Desensitisers

Desensitiser

Contains Potassium Oxalate

- Blocks dentinal tubules.
- For immediate relief from sensitivity.
- Use after vital crown cutting.



Available Packing • **Regular Pack** : 3g syringe & 5 application needles
 • **Twin Pack** : 2 x 3g syringe & 10 application needles

Restorite Desensitiser



Liquid Solution for Hyper Sensitivity

- To treat pain in exposed cervical areas not indicated for restoration
- To alleviate or prevent dentinal sensitivity after preparation of teeth to receive direct or indirect restorations.
- For total etch system Restorite Desensitiser can be used as a re-wetting agent.

Available Packing : 5ml bottle

Prophylaxis & Polishing System

Aster Compo



Alumina Based Polishing Paste for Polishing Dental Composites.

- Non-splattering.
- Thixotropic nature.
- Mint fragrance.
- Fine grit.

Available Packing : 4g Syringe

REF : ATC001

Aster Prophy Prophylaxis Paste



- Dental Prophylaxis Paste for polishing & removal of plaque & stains.
- Thixotropic nature.
- Non-splattering effect.
- Pumice based.
- Fluoride free.



Aster Prophy (Pack of 4)



New Flavours :



Available Packing : 45g tube / Assortment pack of 4 / 200g Jar

REF : ATP001
 REF : ATP002

Miscellaneous

BiofilmX



Dental Chair Unit Waterline Biofilm Reducer

- Ready-to-use solution
- Economical method to maintain dental chair unit tubing hygiene
- 1 bottle offers 4 uses of 125ml each
- For new clinical setup only 125ml once a week is sufficient to obtain optimum results
- For old clinical setup 125ml needed for each day for the first 5 days and there after 125ml once a week is sufficient to obtain optimum results

Available Packing : 500ml Bottle

REF : BFL001

Debubblizer



Reduce Surface Tension Non Foaming Solution

- Spray to be Used Before Pouring Impressions or Wax Patterns
- Helps to avoid bubbles while pouring models.
- Fast drying surface tension reducer & cleanser.

Available Packing : 100ml Spray Bottle / 500ml Refill Jar

REF : DBS001

NT Fog



Spray for Anti-Fog

- Effectively used on MM tops to make it fog free.

Available Packing : 100ml Spray Bottle

REF : NFG001

Caries Indicator



- Dye which shall stain carious dentin to bright red. This would enable dentist to distinguish carious lesion.
- Available colors : Red & Blue.

Available Packing : 15ml bottle

Dent Dry



Dentin Drying & Cleansing Agent

- For having moist dentin surface
- Use before bonding of dentin
- Use before root canal sealant

Available Packing : 30ml bottle


ProsthoClean 2026 Launch



Cleaning Bonding Surface of Prosthesis

- Extra Oral Use Only
- Cleans Bonding Surfaces of All Types of Restorations After Try-In
- Can Be Used with Zirconia, Ceramic, Metal, Composite, etc.
- Optimizes Adhesion with Cements

Available Packing : 5ml bottle



Prime Dental Products Pvt Ltd

Plot No. A-148, Road No. 24, MIDC, Wagle Industrial Estate, Thane (West)-400604, Maharashtra, India.

+91 9272235922 info@prime-dental.in www.prime-dental.in



WhatsApp

